



AccuroFit Studio Companion

GETTING STARTED



SCAN TO DOWNLOAD

STEP 1

Scan QR Code

Scan QR code with your phone camera.



STEP 2

Download AccuroFit Studio Companion app

Download the AccuroFit Studio Companion app to your Apple Watch®



Create your account


Tap 'Sign in' to create your account.



STEP 4

Setup your Sensor ID

Enter your existing Sensor ID or generate a new one.



STEP 5

Connecting to the Hub

The Watch app will automatically connect to your Accurofit hub.



STEP 6

Start your workout

Tap to start your workout. You should see your heart rate on the AccuroFit system.

TIPS & TRICKS

USE ONLY STUDIO HUB APP

Avoid running Apple's workout app at the same time as the *AcurroFit Studio Companion app*, it will not work.



SAVING TO APPLE HEALTH

The *AcurroFit Studio Companion app* can save all of the same information to Health as using Apple's workout app with the added ability to connect to your studio's group fitness system.

Be sure the 'Save to Health' toggle is enabled in the app's settings.

FILLING YOUR RINGS

When enabled, the *AcurroFit Studio Companion app* will save all of your workout history and data to Apple Health including contributing to filling your rings!

Open the Activity app on the iPhone to view your complete workout history and data for your workouts.



SELECTABLE WORKOUT TYPES

The *AcurroFit Studio Companion app* can save a variety of workout types when saving to Health.

Open the app's settings to select the appropriate workout type for your studio's workout.

